

# System Serial Number Locations for PPMS, DynaCool and VersaLab systems

## 1) PPMS, PPMS EverCool, and PPMS EverCool-II:

(Example S/N PPMS-603)



(wooden cabinet models only)  
Front of PPMS Cabinet.  
Open door at bottom.  
Lower left side.



(all cabinet models)  
Back of PPMS Cabinet.  
Lower right corner.

## 2) DynaCool:



**Back of DynaCool cabinet.**

**Upper left, at the top of the electronics rack.**

## 3) VersaLab:



**Back of VersaLab unit.**

**Lower right, near vent gas port.**