




Service Note 1070-211

How to Locate the Serial Number on Quantum Design Systems

- I. Magnetic Property Measurement Systems (MPMS)
- II. Physical Property Measurement Systems (PPMS)
- III. MPMS Squid VSM (MPMS3)
- IV. VersaLab
- V. PPMS DynaCool
- VI. Advanced Technology Liquefier (ATL)
- VII. Quantum Design Compressor

I. Magnetic Property Measurement System (MPMS)

1. The serial number sticker is located on the lower right rear of the cabinet.



II. Physical Property Measurement System (PPMS)

Wooden Cabinet

1. Open the front door at the bottom of the cabinet. The serial number is located on the lower left side.



2. On the lower right rear of the cabinet is a second serial number sticker.




Metal Cabinet

1. The serial number sticker is located on the lower right rear of the cabinet.



III. MPMS Squid VSM (MPMS3)

1. The serial number sticker is located on the right rear of the cabinet.




IV. VersaLab

1. The serial number sticker is located on the lower right rear of the cabinet.




V. PPMS DynaCool

1. The serial number sticker is located on the top left rear of the cabinet.



VI. Advanced Technology Liquefier (ATL)

1. The serial number sticker is located on the top left of the ATL Unit.



VIII. Quantum Design Compressor

1. The serial number is located on the lower left rear of the cabinet.

