

PPMS Family Configurations

Legend:		QD Measurement Option							Other Inserts		
		magnetic			thermal		Electrical				
Available		VSM	ACMS	Torque	Heat Cap	TTO	ETO	ACT	Resistivity		
Hydrostatic pressure cell											
Magneto-Optic Option											
Not available											
platform	insert										
PPMS	platform	Available	Available	Not available	Available	Available	Available	Available	Available	Available	SPM ADR FMR
	Rotator	Not available	Not available	Available	Not available	Not available	Available	Available	Available	Available	
	³He	Not available	Not available	Available	Available	Not available	Available	Available	Available	Available	
	Dilution	Not available	Available	Not available	Available	Not available	Available	Not available	Not available	Not available	
PPMS EC-II	platform	Available	Available	Not available	Available	Available	Available	Available	Available	Available	ADR FMR
	Rotator	Not available	Not available	Available	Not available	Not available	Available	Available	Available	Available	
	³He	Not available	Not available	Available	Available	Not available	Available	Available	Available	Available	
DynaCool	platform	Available	Available	Not available	Available	Available	Available	Available	Available	Available	ADR FMR
	Rotator	Not available	Not available	Available	Not available	Not available	Available	Not available	Available	Available	
	³He	Not available	Not available	Available	Available	Not available	Available	Available	Available	Available	
	Dilution	Not available	Available	Not available	Available	Not available	Available	Not available	Not available	Not available	
VersaLab	platform	Available	Available	Not available	Available	Available	Available	Available	Available	Available	FMR
	Rotator	Not available	Not available	Available	Not available	Not available	Available	Not available	Available	Available	

Please consult sales representative for option availability